



SIGVARIS
GROUP

ASSORTMENT
PHLEBOLOGY 2023



COLORFUL EVERYDAY – TREND COLORS 2023/2024

What a menu! Serving you our new TREND COLORS: salted caramel, blueberry dream, fresh basil and hot chili.

What a feast. Not only for the eyes, but also for the legs. The fresh colors from our popular Style SEMITRANSSPARENT line make legs look good and give an elegant, fashionable and exciting appearance. Colors to fall in love with.



COMPLETE PRODUCT RANGE ACCORDING TO COMPRESSION LEVEL

Product	Model	Rim	CCL				Toe		SVGT*	Grip Top		Page	
			1	2	3	4				KGT*	MKGT*		
sigvaris													
Essential	COMFORTABLE	A-D	Soft top	1	2			•	•			13	
			Grip top	1	2			•	•		•		
		A-G	Soft top	1	2			•	•				
			Grip top	1	2			•	•	•	•		
			Waist attachment		2			•	•				
		A-T	Soft top	1	2			•	•				
	A-T Maternity	Soft top	1	2			•	•					
	A-T Bodyform	Soft top		2			•	•					
	THERMO-REGULATING	A-D	Soft top	1	2	3		•	•			15	
			Grip top	1	2			•	•		•		
		A-G	Soft top	1	2			•	•				
			Grip top	1	2	3		•	•	•	•		
			Waist attachment	1	2	3		•	•	CCL 1 / 2			CCL 1 / 2
		A-T	Soft top	1	2	3		•	•				
	A-T Maternity	Soft top	1	2			•	•					
	CLASSICAL	A-D	Soft top	1	2			•	•			17	
			Grip top	1	2			•	•		•		
A-F		Grip top		2			•	•	•	•			
A-G		Grip top	1	2			•	•	•	•			
		Waist attachment	1	2			•	•					
A-T		Soft top	1	2			•	•					
A-T Maternity	Soft top		2			•	•						
Style	SEMI-TRANSPARENT	A-D	Soft top	1	2			•	•			21	
			Grip top	1	2			•	•		•		
		A-G	Grip top	1	2			•	•	•	•		
		A-T	Soft top	1	2			•	•				
A-T Maternity	Soft top		2			•	•						
Active Specialities	MASCULINE	A-D	Soft top		2				•			25	
			Grip top		2								
	TRADITIONAL	A-D	Soft top		2	3	4		•			29	
			Grip top		2	3			•		•		
		A-F	Soft top		2	3	4		•				
			Grip top		2	3			•		•		
		A-G	Soft top		2	3	4		•				
			Grip top		2	3			•		•		
			Waist attachment		2	3			•				
		A-G Post-OP	Waist attachment		2	3			•				
	A-T	Soft top		2				•					
	A-T Maternity	Soft top		2				•					
	A-T For men	Soft top		2				•					
	ULCER X	A-D	Soft top								Understocking closed, TRADITIONAL open	31	
	DIABETIC COMPRESSION SOCKS	A-D	Soft top							•		33	
	THROMBO / THROMBO X	A-D	Soft top									35	
			Grip top								•		
		Waist attachment											
ARM SLEEVE ADVANCE	A-G with mitten	Soft top									37		
		Grip top							•				
	C-G without mitten	Soft top											
Grip top								•					
ARM SLEEVE TRADITIONAL	A-C only mitten	-		2									
	C-H without mitten	Support straps		2									
	C-G without mitten	Grip top		2						•			
	A-H with mitten	Support straps		2									
	A-G with mitten	Grip top		2						•			
DYNAVEN													
Dynaven	CLASSIC	A-D	Soft top	1	2			•	•			47	
			Grip top	1	2			•	•		•		
		A-G	Grip top	1	2			•	•		•		
		A-T	Soft top	1	2			•	•				
		A-T Maternity	Soft top	1	2			•	•				

* SVGT = Sensinova grip top; KGT = Quadrinova knobbed grip top (Essential COMFORTABLE, Essential THERMOREGULATING, Essential CLASSICAL, Style SEMITRANSPARENT, Dynaven CLASSIC) / Knobbed grip top (Specialities TRADITIONAL, THROMBO / THROMBO X and Arm Sleeve TRADITIONAL); MKGT = Micro knobbed grip top












COMPLETE PRODUCT RANGE ACCORDING TO INDICATIONS / NEEDS

Indication / Need	Segment / Product name	Page
Calf stocking with more grip on the leg: Calf with grip top	Essential COMFORTABLE	13
	Essential THERMOREGULATING	15
	Essential CLASSICAL	17
	Style SEMITRANSSPARENT	21
	Specialities TRADITIONAL	29
	Dynaven CLASSIC	47
Voluminous feet or feet with overlength	Essential THERMOREGULATING (maxifoot: foot with overlength)	15
	Active MASCULINE (voluminous foot)	25
Shorter thigh length: Half-thigh	Essential CLASSICAL	17
	Specialities TRADITIONAL	29
Thigh stocking for sensitive / allergy-prone skin: Thigh with soft top	Essential COMFORTABLE	13
	Essential THERMOREGULATING	15
	Specialities TRADITIONAL	29
Thigh stocking variant for easier donning and doffing: Thigh with waist attachment	Essential COMFORTABLE	13
	Essential THERMOREGULATING	15
	Essential CLASSICAL	17
	Specialities TRADITIONAL	29
Short lower leg or leg length	Essential THERMOREGULATING (A-D Calf and A-T Pantyhose)	15
Beginning lipedema in the hip area / soft connective tissue: Pantyhose Bodyform	Essential COMFORTABLE	13
Ulcerated leg	Specialities ULCER X	31
Diabetic foot	Specialities DIABETIC COMPRESSION SOCKS	33
Postoperative care: Post-OP	Specialities TRADITIONAL	29
Distinctive edema / highly advanced CVI: Compression class 3 and/or 4	Essential THERMOREGULATING	15
	Specialities TRADITIONAL	29
Thrombosis prophylaxis	Specialities THROMBO / THROMBO X	35
Prevention / treatment and care of lymphedema of the upper extremities	Specialities ARM SLEEVE ADVANCE	37
	Specialities ARM SLEEVE TRADITIONAL	37
Sensitive skin	Essential THERMOREGULATING	15
Elegant and semitransparent stocking	Style SEMITRANSSPARENT	21
Thermoregulating effect	Essential THERMOREGULATING	15

Explanation on special models / special requirements on **p.50**

COMPLETE PRODUCT RANGE ACCORDING TO MODELS

Product		 A-D Calf	 A-F Half-thigh	 A-G Thigh	 A-G with waist attachment	 A-G Post-OP	 A-T Pantyhose	 A-T Maternity	 A-T Bodyform	 A-T For men	Page
Essential	COMFORTABLE	•		•	•		•	•	•		13
	THERMOREGULATING	•		•	•		•	•			15
	CLASSICAL	•	•	•	•		•	•			17
Style	SEMITRANSSPARENT	•		•			•	•			21
Active	MASCULINE	•									25
Specialities	TRADITIONAL	•	•	•		•	•	•		•	29
	ULCER X	•									31
	DIABETIC COMPRESSION SOCKS	•									33
	THROMBO / THROMBO X	•		•	•						35
	ARM SLEEVE ADVANCE	available as A-G with mitten and C-G without mitten									37
	ARM SLEEVE TRADITIONAL	available as A-C only mitten, C-H and C-G without mitten as well as A-H and A-G with mitten									37
Dynaven	CLASSIC	•		•			•	•			47

MATERIAL COMPOSITION

Product	Model / Material composition	
Essential	COMFORTABLE	all models: 65% polyamide, 35% elastane
	THERMOREGULATING	all models: 52% polyamide, 34% elastane, 14% cotton
	CLASSICAL	A-D, A-F, A-G: 68% polyamide, 32% Elasthan; A-T, A-T Maternity: 76% polyamide, 24% elastane
Style	SEMITRANSSPARENT	A-D: 61% polyamide, 39% elastane; A-G: 63% polyamide, 37% elastane; A-T, A-T Maternity: 59% polyamide, 41% elastane
Active	MASCULINE	all models: 81% polyamide, 19% elastane
Specialities	TRADITIONAL	all models: 70% polyamide, 30% elastodiene
	ULCER X	Cover stocking: see Specialities TRADITIONAL; Understocking: 62% polyamide, 26% elastane, 12% cotton
	DIABETIC COMPRESSION SOCKS	all models: 65% polyester, 19% polyamide, 11% cotton, 5% elastane
	THROMBO / THROMBO X	THROMBO: 74% polyamide, 26% elastane; THROMBO X: 81% polyamide, 17% elastane, 2% polypropylene
	ARM SLEEVE ADVANCE	all models: 47% polyamide, 33% Modal (Micro-Modal), 20% elastane
	ARM SLEEVE TRADITIONAL	all models: 70% polyamide, 30% elastodiene
Dynaven	CLASSIC	A-D, A-G: 68% polyamide, 32% Elasthan; A-T, A-T Maternity: 76% polyamide, 24% elastane

TABLE OF CONTENT

About SIGVARIS GROUP	2
Why SIGVARIS GROUP	3
Our customer groups	4
Services and labels	5
Colors	6
Grip tops	7

sigvaris 8

Essential	11
COMFORTABLE	13
THERMOREGULATING	15
CLASSICAL	17

Style	19
SEMITRSPARENT	21

Active	23
MASCULINE	25

Specialities	27
TRADITIONAL	29
ULCER X	31
DIABETIC COMPRESSION SOCKS	33
THROMBO / THROMBO X	35
ARM SLEEVES ADVANCE & TRADITIONAL	37

Custom made and alterations	38
------------------------------------	-----------

Accessories	40
DONNING AND DOFFING DEVICES	40
CARE AND ADDITIONAL PRODUCTS	42

DYNAVEN 45

Dynaven CLASSIC	47
-----------------	----

Information	48
Taking measurements correctly	48
Tips and tricks	49
Special models/special requirements	50
Compression values, indications and contraindications	52


Sizing systems	54
-----------------------	-----------

ABOUT SIGVARIS GROUP

We want to develop the best compression solutions for people. Each and every individual should feel good. Every day. We have been focusing on this for over 150 years.

Our portfolio is made up of various segments that address various needs and indications. Each segment is designed to not only let you find a suitable, but also the best compression solution for your customer. On the one hand, the focus is on the effect and medical aspects – on the other, the products should make the customer feel good and help to lead a life without restrictions.

Quality made in Switzerland is important to us. Therefore we are extremely proud of our production facility in St. Gallen.



WHY SIGVARIS GROUP

Our compression products fit into every lifestyle, as we focus on people. This is also reflected in our product range. We have understood that people simply want to be themselves. No matter where and when.

You know your customers. And your customers trust you to help them feel good. The newly structured range makes it easier to find the best solution. We have arranged our proven products into segments to help you to find the right product with ease.

WE ARE DEDICATED TO HELPING PEOPLE FEEL THEIR BEST. EVERY DAY.

OUR SEGMENTS

Essential

These solutions are designed for life situations where customers are looking for something uncomplicated. Comfort and reliability are the top priorities. For everyday suitability.

Style

The fashion aspect is the focus here. Style is guaranteed to offer the right stocking for every occasion – be it everyday office life, a relaxing brunch or a special event.

Active

For a lifestyle with lots of physical activity – whether in your leisure time or at work. Active offers the perfect solution for active days.

Specialities

When medical aspects dominate and standard compression solutions are not sufficient: these products are intended for very specific medical indications or needs.

Dynaven

The functional and effective stocking for pragmatic users. No more and no less.



OUR CUSTOMER GROUPS

Every person is unique and their requirements for compression solutions can change depending on their life situation. For example, a product from the Essential segment may prove to be the right one for during the day, and for an evening out, Style may be the one that fits best. Or someone may decide to use Active during their daily work, but chooses Essential for relaxing hours at home.

The structure of the range is tailored to the needs of the people who use our products. In a comprehensive market research study we interviewed them and then consistently oriented ourselves according to the results.



The Essential clientele
For them, the main thing is that their product is uncomplicated and convenient. Their segment is **Essential**.



The Style clientele
For them, fashion aspects are important. Attractive, stylish stockings are important. They can find them in the **Style** segment.



The Active clientele
Their dynamic lifestyle should also be reflected in the compression solutions. The products should be robust and provide active support, just like the products in the **Active** segment.



The Dynaven clientele
They attach importance to functionality and effectiveness. They can find their stockings under **Dynaven**.

SERVICES AND LABELS

Offering you the highest service and product quality – that is important to us. We stand for compression therapy and are pleased to provide support with our product knowledge and medical expertise.

In addition to our internal quality assurance, which meticulously tests our products, various quality seals also oblige us to maintain high quality.



We are involved in our products and compression therapy on a daily basis. We pass this knowledge on through training, our field force and customer service. We are constantly working on optimizing our services such as consulting services and ordering processes. It is also important to us that we make compression therapy known to a broader public.



We fulfill all quality and test regulations of the RAL Deutsches Institut für Gütesicherung und Kennzeichnung e.V. (German Institute for Quality Assurance and Labeling) for medical compression stockings.




Our medical compression stockings have been tested for human ecology and comply with the internationally binding Eco-Tex standard 100.



Our medical compression stockings meet the general safety and performance requirements of the relevant EU regulations.

COLORS

Today compression stockings are available in a wide range of colors. In addition to bold fashionable colors, numerous standard colors – such as various skin tones and black – are also available. The stockings in the Style segment in particular are characterized by a large selection of colors. The Style SEMITRANSSPARENT range features four new trend colors every two years.




GRIP TOPS

Depending on the product, the wearers of compression stockings can choose from different grip top options. All grip tops impress with their skin compatibility and high wearing comfort.

Essential COMFORTABLE
Essential THERMOREGULATING
Essential CLASSICAL



Sensinova grip top SVGT
neutral design; invisible



Quadrinova knobbed grip top KGT
neutral design; pleasant against the skin

Essential THERMOREGULATING



Micro knobbed grip top MKGT
neutral design; for sensitive skin

Style SEMITRANSSPARENT




Sensinova grip top SVGT
feminine design; invisible



Quadrinova knobbed grip top KGT
feminine design; pleasant against the skin

Specialities TRADITIONAL



Knobbed grip top KGT
neutral design

Dynaven CLASSIC



Quadrinova knobbed grip top KGT
neutral design; pleasant against the skin

Live your life. Be yourself.

Sigvaris compression stockings meet the requirements of wearers who expect more than mere functionality and effectiveness from their chosen product. In the Essential, Style, Active and Specialities segments, discerning users will find products that fit their specific needs and personal preferences and are suitable for every moment of their day. All products convince through individual features and high quality workmanship.

Welcome to the world of Sigvaris.





Essential

Comfort and reliability are the top priorities in the Essential Line solutions. They are the flexible everyday companions, uncomplicated and yet visually appealing. They fit the female leg just as well as the male leg. Each model is designed for regular, uncomplicated and reliable use in the pursuit of health.

Essential COMFORTABLE p. 13

Essential THERMOREGULATING p. 15

Essential CLASSICAL p. 17



"I choose Essential for my daily work routines. Simple donning and doffing and uncomplicated handling are my top priorities."

"I prefer reliable and comfortable compression solutions no matter what I do. That is why I choose Essential."





Essential COMFORTABLE

For her and him

Essential COMFORTABLE is the compression stocking that fits like a glove and meets the highest standards of comfort and care. It is comfortably soft, pleasant to wear and adapts perfectly to one's shape.

COMFORTABLE		A-D Calf		A-G Thigh		A-T Pantyhose	A-T Maternity	A-T Bodyform	
Matchcode ECF		Soft top	Grip top	Soft top	Grip top	Waist attachment	Soft top	Soft top	Soft top
CCL	1	•	•	•	•		•	•	
	2	•	•	•	•	•	•	•	•
Toe		•	•	•	•	•	•	•	•
		•	•	•	•	•	•	•	•
Grip top	Sensinova grip top				•				
	Quadrinova knobbed grip top		•		•				
Colors	black	•	•	•	•	•	•	•	•
	savanna	•	•	•	•	•	•	•	•
	caramel	•	•	•	•	•	•	•	•
	skin	•			•		•	•	
	grey*	•			•		•	•	

* No stock item, produced on order

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models/needs on **p.50**



perfect fit
Comfortable: perfect fit for her and him



opaque
Attractive: non-transparent, conceals venous problem zones perfectly



easy to use
Convenient: easy to put on and take off










Essential THERMOREGULATING

For her and him

The fine cotton layer of Essential THERMOREGULATING lies directly on the skin and has a thermoregulating effect. The cotton absorbs excess moisture directly from the skin and stores it inside the fibers. This prevents developing a moisture film on the skin which promotes the undesirable growth of bacteria. This is beneficial for people with sensitive skin.

THERMOREGULATING		A-D Calf		A-G Thigh			A-T Pantyhose	A-T Maternity
Matchcode ETH		Soft top ^{1,2}	Grip top	Soft top	Grip top	Waist attachment	Soft top ³	Soft top
CCL	1	●	●	●	●	●	●	●
	2	●	●	●	●	●	●	●
	3	●			●	●	●	
Toe		●	●	●	●	●	●	●
		●	●	●	●	●	●	●
Grip top	Sensinova grip top				● CCL 1/2			
	Quadrinova knobbed grip top		●		●			
	Micro knobbed grip top				● CCL 1/2			
Colors	 black	●	●	●	●		●	●
	 nature	●	●	●	●	●	●	●
	 marine*	available as A-D with soft top, CCL 2, closed toe, normal and long						

* No stock item, produced on order

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models/needs on **p. 50**

¹ CCL 2 with closed toe is also available as **maxifoot** in black or nature, without XXS and XXS Plus

² CCL 2 with open and closed toe, is also available as short length in black or nature

³ CCL 2 and 3 with open and closed toe, is also available as short length in nature



sensitive skin

Soft: is also suitable for sensitive legs due to 14 % cotton content and soft top



thermoregulating

Cooling: evaporation of moisture stored in the cotton, resulting in a cooling effect

Warming: the air in the spaces between the knitted fabric and the fibers forms an insulating layer at cold temperatures



opaque

Attractive: discreet, non-transparent ribbed structure, conceals venous problem zones perfectly





Essential CLASSICAL

For her and him

With its non-transparent appearance, the Essential CLASSICAL is the discreet companion for everyday situations. The integrated microfiber ensures pleasant wearing comfort and a discreet appearance.

CLASSICAL		A-D Calf		A-F Half-thigh	A-G Thigh		A-T Panythose	A-T Maternity
Matchcode ECL		Soft top	Grip top	Grip top	Grip top	Waist attachmet	Soft top	Soft top
CCL	1	●	●		●	●	●	
	2	●	●	●	●	●	●	●
Toe		●	●	●	●	●	●	●
		●	●	●	●	●	●	●
Grip top	Sensinova grip top			●	●			
	Quadrinova knobbed grip top		●	●	●			
Colors	black	●	●	●	●	●	●	●
	savanna	●	●	●	●	●	●	●
	caramel	●	●	●	●	●	●	●
	crispa	●	●	●	●	●	●	●

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models/needs on **p.50**



easy to use

Convenient: easy to put on and take off



opaque

Attractive: non-transparent, conceals venous problem zones perfectly



pleasant wearing comfort

Comfortable: high wearing comfort



Style

Sigvaris Style not only focuses on the medical aspect but also on the fashion aspect. Here wearers are able to express the different facets of their personality and underline their individuality. Style offers the right stocking for everyday life, as well as for special occasions. The models are designed for an elegant, self-confident and attractive appearance. The stockings convince with sophisticated product details and go just as well with evening wear as with business outfits.

Style SEMITRANSSPARENT p. 21



"I love fashion and colors. Attractive product details are also important to me, be it in the office, at home or at a party. I can find my compression stocking in the Style segment."

















Style SEMITRANSSPARENT

For her

Style SEMITRANSSPARENT fascinates with stylish, elegant charisma, innovative details and maximum wearing comfort. It meets the highest demands of style and fashion and convinces with its semitransparent appearance. The compression stocking was developed specifically for female needs and impresses with its high wearing comfort and that extra soft feel.

SEMITRANSSPARENT		 A-D Calf		 A-G Thigh		 A-T Pantyhose		 A-T Maternity	
Matchcode SSE		Soft top	Grip top	Grip top	Soft top	Soft top			
CCL	1	•	•	•	•	•			
	2	•	•	•	•	•			
Toe		•	•	•	•	•			
		•	•	•	•	•			
Grip top	Sensinova grip top			•					
	Quadrinova knobbed grip top		•	•					
Colors	 black  savanna  caramel	•	•	•	•	•			
	 skin	•		•	•	•			
No stock item, produced on order:	 grey  marine								
	 salted caramel  blueberry dream	•		•	•	•			
	 fresh basil  hot chili								

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on **p.50**



semitransparent

Fine: elegant, semitransparent appearance



elegant

Stylish: feminine esthetics; also available in four trendy colors



perfect shaping

Attractive: shapes leg and buttocks optimally





"My compression stocking has to keep up with my active life and support me and my legs. That is why I chose Active."



Active

Sigvaris Active stockings are suitable for people who are dependent on special support in their dynamic lives. These solutions ensure optimum support, even in active situations. The product properties ensure that the wearer feels confident and safe at any time of the day or year.






Active MASCULINE p. 25



Active MASCULINE

For him

Active MASCULINE is the perfect solution for all men who lead active lives. Due to moisture-regulating yarns and the integrated odor-inhibiting zinc finish, feet stay fresh even in critical situations. Active MASCULINE was developed especially for men's feet, which guarantees an optimal fit.

MASCULINE		
Matchcode	AMA	
CCL	2	●
Toe		●
Colors	 black	●
	 carbon	●
	 cyanide	●

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on **p.50**



cushioning sole

Comfortable: with full-surface sole padding



reinforced heel

Strong: reinforced heel for more safety in active everyday life



seamless toe section

Pleasant to wear: seamless toe for fewer pressure points



Specialities

Sometimes the medical aspect is of primary concern and the needs of the stocking wearer cannot be met with standard compression solutions. These are the situations where products from the Specialities segment come into their own. They are characterized in that each stocking has been developed for a very specific medical indication or special need.

Specialities TRADITIONAL	p. 29
Specialities ULCER X	p. 31
Specialities DIABETIC COMPRESSION SOCKS	p. 33
Specialities THROMBO / THROMBO X	p. 35
Specialities ARM SLEEVES ADVANCE & TRADITIONAL	p. 37



Specialities TRADITIONAL

For her and him

Whenever high medical effectiveness is required, TRADITIONAL is the right stocking. The unique stress-strain behavior of the yarn composition ensures high pressure stability throughout the entire day and effective treatment even in advanced venous diseases. TRADITIONAL is easy to don even at high compression levels. It is robust, durable and breathable.

TRADITIONAL		A-D Calf		A-F Half-thigh		A-G Thigh		A-G Post-OP*	A-T Pantyhose	A-T Maternity	A-T For men	
Matchcode TRA		Soft top	Grip top	Soft top	Grip top	Soft top	Grip top	Waist attachment	Waist attachment	Soft top	Soft top	Soft top
CCL	2	●	●	●	●	●	●	●	●	●	●	●
	3	●	●	●	●	●	●	●	●			
	4	●		●		●						
Toe		●	●	●	●	●	●	●	●	●	●	●
Grip top	Knobbed grip top		●		●		●					
Colors	beige	●	●	●	●	●	●	●	●	●	●	●

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models/needs on **p.50**

* Kit comprises: 2x understocking (without compression), 1x TRADITIONAL

Available sizes on **p.57**

Reliable: high safety and medical efficacy for successful therapy

Robust: maximum compression stability throughout the day due to refined natural rubber, excellent force-elongation properties



Effective: recommended by doctors worldwide





Specialities ULCER X

For her and him

ULCER X is a compression stocking system for the treatment of ulcerated legs (ulcus cruris venosum) consisting of a liner (15 - 20 mmHg; 2.0 - 2.7 kPa) and a cover stocking TRADITIONAL (CCL 2). Compared to bandages, ULCER X is easy to put on, usually without external assistance, and comfortable to wear. The compression stocking system allows controlled reproducible compression.

ULCER X		
Matchcode ULX		
CCL	~45mmHg (~6.0 kPa)	●
Toe		● (TRADITIONAL)
		● (Understocking)
Colors	● beige	●

Studies attribute ULCER X various advantages – for more information, please ask your customer service or field force representative for the ULCER X folder.

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on **p.50**

Available sizes on **p.57**

Protection: due to the patented liner, the wound dressing also remains fixed in place at night and the wound is not affected adversely when used correctly – neither when donning nor doffing

Effective: better healing rates in small and medium-sized ulcers compared to bandages



ULCER X Kit comprises: 2x understocking, 1x TRADITIONAL
Additional understockings are available separately (4 pieces)
Notes on size of understockings on **p. 57**







Specialities

DIABETIC COMPRESSION SOCKS

For her and him

DIABETIC COMPRESSION SOCKS are suitable for people with diabetes with an indication for compression due to diabetes or venous lower leg edema or symptoms of vein dysfunction. The sock responds to the special needs of diabetic feet. It features light compression pressure (18 - 25 mmHg; 2.4 - 3.3 kPa), a compression-free foot section, soft knitted fabric and special padding on the sole of the foot and in the toe area. A seamless toe has been specifically developed to protect the toes.

DIABETIC COMPRESSION SOCKS		
Matchcode DCS		A-D Calf
		Soft top
CCL	18-25 mmHg (2.4-3.3 kPa)	<input checked="" type="checkbox"/>
Toe		<input checked="" type="checkbox"/>
Colors	<input checked="" type="checkbox"/> white	<input checked="" type="checkbox"/>

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on **p. 50**

It is recommended to use DIABETIC COMPRESSION SOCKS following a medical examination. For detailed information on use as well as on stockings, please refer to the explanations in the product communication or the package insert.

A special men's and women's variant guarantees the best possible fit in the foot area. The two variants differ in the size of the foot section.

Contraindications for DIABETIC COMPRESSION SOCKS WOMEN / MEN

Medical compression stockings must not be worn in the following cases: Severe arterial circulatory disorder (if any of these parameters apply: ABI <0.5, ankle artery pressure <60 mmHg, toe pressure <30 mmHg, or TcPO2 <20 mmHg dorsum of foot); decompensated heart failure; advanced peripheral neuropathy; intolerance to compression garment materials; pronounced, weeping dermatoses; pronounced edema in the foot area and primary chronic polyarthritis. Relative contraindications, whereby benefits and risks must be weighed up, in particular: arterial circulatory disorders (ABI between 0.5 and 0.8).

Skin-caring: water-absorbing and water-repellent yarn components provide a good moisture climate in the sock

Safety: white color gives a clear view of any weeping wounds

Hygiene: washable at 95 °C in the laundry net



Specialities THROMBO / THROMBO X

For her and him

The THROMBO / THROMBO X is used for thrombosis prophylaxis in recumbent patients. It is suitable following surgical interventions (especially abdominal, lower abdominal and orthopedic surgery), after fractures or injuries to the legs and in cases of restricted mobility.

THROMBO / THROMBO X		A-D Calf		A-G Thigh	
Matchcode THR / THX		Soft top		Grip top	
				Waist attachment	
CCL	~18 mmHg (~2.4 kPa)	●	●	●	●
Toe	with inspection hole	●	●	●	●
Grip top	Knobbed grip top		●		
Colors	○ white	●	●	●	●

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on **p.50**

* Size XS not available for THROMBO A-G Thigh with waist attachment

Gentle: degressive pressure curve and compression of ~ 18 mmHg (~ 2.4 kPa) in the ankle area

Caring: skin-friendly and breathable

Practical: inspection hole available on all models






Specialities

ARM SLEEVES ADVANCE & TRADITIONAL

For her and him

Arm compression sleeves are used for the prevention as well as treatment and care of lymphedema of the upper extremities and are an important part of maintenance therapy.

ARM SLEEVE ADVANCE		 A-G with mitten		 C-G without mitten	
Matchcode AAD		Soft top*	Grip top	Soft top*	Grip top
CCL	14–18 mmHg (1.8–2.4 kPa)	●	●	●	●
	20–25 mmHg (2.7–3.3 kPa)	●	●	●	●
Grip top	Sensinova grip top		●		●
Colors	 light beige	●	●	●	●







*no Plus sizes available

Skin-caring: the polyamide microfiber with silver-based additive offers a soft texture.

In addition, it inhibits the growth of bacteria and prevents odors and skin irritations

Wear-friendly: the thermoregulating modal fiber absorbs moisture particularly well and is therefore extremely comfortable to wear



ARM SLEEVE TRADITIONAL		 A-C Only mitten	 C-H without mitten with support straps	 C-G without mitten with grip top	 A-H with mitten with support straps	 A-G with mitten with grip top
Matchcode ATR						
CCL	2	●	●	●	●	●
Grip top	Knobbed grip top			●		●
Colors	 beige	●	●	●	●	●

Skin-friendly: high breathability due to large-mesh knit

Pressure stability: high recovery force due to natural rubber

Wear-resistant: up to 6 times better* abrasion results (* = less wear) in direct comparison with jeans (proven in wear tests)



Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on p.50

CUSTOM MADE AND ALTERATIONS

Our product segments cover many customer needs as standard. However it may happen that a standard product cannot meet the special requirements of a customer. In this case, you can resort to our made to measure range and our alteration service.

CUSTOM MADE ASSORTMENT

- Perfect fit– precisely tailored to the personal leg measurements
- The sign of bespoke work – the customer's name can be found on the stocking label
- 100 % satisfaction guaranteed

Please ask for our order form for custom-made stockings. Our customer service will be pleased to answer any questions you may have about our range of made to measure products.

Model	Top	Essential COMFORTABLE	Essential THERMO-REGULATING	Essential CLASSICAL	Style SEMITRANS-PARENT	Specialities TRADITIONAL
A-D Calf	Soft top	●	●	●	●	●
	KGT (narrow / broad)	●	●	●	●	●
	SVGT	●	●	●	●	●
A-F Half-thigh	Soft top	●	●	●		● ³
	KGT broad	●	●	●	● ²	● ³
	SVGT	●	●	●	● ²	● ³
A-G Thigh	Soft top	●	●	●		● ³
	KGT broad	●	●	●	●	● ³
	MKGT		● ¹			
	SVGT	●	●	●	●	● ³
A-T Pantyhose	Soft top		●	●	●	● ³
		●	●	●	●	● ³
A-T For men	Soft top		●	●		● ³
		●	●	●		● ³
A-T Maternity	Soft top		●	●	●	● ³
		●	●	●		● ³
A-T Bodyform	For women / men		●	●		●
	With fly		●			●
A-T one leg pantyhose	For women / men		●	●		● ³
		●	●	●		● ³
	CCL	1 2	1 2 3	1 2	1 2	2 3 4
	Colors	savanna, skin, grey ⁴ , caramel, black	nature, black	caramel, crista, savanna, black	savanna, skin, grey ⁴ , marine ⁴ , caramel, black, trend colors ⁴	beige

¹ only available in CCL1 and CCL2 ² only available in CCL2
³ compression up to the cG level (without soft-part) ⁴ No stock item, produced on order

Available grip tops:
 KGT: Knobbed grip top (KGT, Quadrinova): 2.5 cm (narrow), 5.0 cm (broad); MKGT: Micro knobbed grip top, 5.0 cm;
 SVGT: Sensinova grip top feminine design, 5.5 cm; Sensinova grip top neutral design, 5.0 cm;
 See overview on grip tops on **p. 7**

All products are available with an open or closed toe.

ALTERATIONS

Our alteration service offers numerous modification options for standard as well as customized products. Special wishes can already be considered with the order.

This only extends delivery times by 3 working days.

Model	Sew on the grip top	Shorten point of the foot [*]	Shortening of the pantyhose part	Including a fly	Cusset change	Stoma opening	Opening in crotch	Shortening the length**
A-D Calf	●	●						●
A-F Half-thigh	●	●						●
A-G Thigh	●	●						●
A-G Thigh with waist attachment		●						●
A-T Pantyhose		●	●	●	●	●	●	
A-T For men with fly (left, right or horizontal)		●	●		●	●	●	
A-T Maternity		●	●		●	●	●	
A-T Bodyform		●	●		●			

* only possible with closed toe
 ** incl. sewing on a new grip top

Available grip tops:
 KGT: Knobbed grip top (KGT, Quadrinova): 2.5 cm (narrow), 5.0 cm (broad);
 SVGT: Sensinova grip top feminine design, 5.5 cm; Sensinova grip top neutral design, 5.0 cm
 See overview on grip tops on **p. 7**

No changes possible for Dynaven CLASSIC as well as Active MASCULINE.

ACCESSORIES

Donning and doffing devices

Many first-time users ask themselves how to put on and take off their compression stockings easily and without great effort. Whether for mobile wearers, or those with limited mobility or for nursing staff - we have suitable solutions for easy donning and doffing of stockings.

FOR EVERY WEARER OF STOCKINGS: Gloves and foot slips

Gloves make it easier to don and doff the stocking and should be handed over with each compression stocking. They help to spread the stocking, protect it and ensure longer durability.

Textile gloves: The latex coating provides a good grip and the knitted fabric is particularly skin-friendly.

Sizes:

- XS
- S
- M
- L
- XL



Rubber gloves: The traditional glove made of rubber for a perfect grip

Sizes:

- S (glove size 6.5–7)
- M (glove size 7.5–8)
- L (glove size 8.5–9)
- XL (glove size 9.5–10)



Foot slips: Facilitates easier donning and doffing of medical compression stockings with an open toe.



FOR MOBILE STOCKING WEARERS OR WHEN TRAVELLING: sim-slide and magnide on/off

These aids are slipped over the foot. The smooth material then helps the stocking to slide across your foot quickly and without any friction.

sim-slide: For compression stockings with open toe.

Sizes:

- S (shoe sizes 34–36)
- M (shoe sizes 37–38.5)
- L (shoe sizes 39–41)
- XL (shoe sizes ≥ 42)



magnide on/off: For compression stockings with closed toe.

Sizes:

- M (shoe sizes ≤ 38)
- L (shoe sizes 39–43)
- XL (shoe size ≥ 44)



FOR STOCKING WEARERS WITH RESTRICTED MOBILITY: MELANY

With MELANY, compression stockings can be donned and doffed quickly, easily and without much effort. Even with limited mobility, this donning and doffing aid can be handled with ease.

Size	cC
XS	25–30 cm
S	30–38 cm
M	38–46 cm
L / L Plus	46–56 cm



cC widest circumference of the calf



ENERGY-SAVING SOLUTION FOR STOCKING WEARERS AND NURSING STAFF: DOFF N' DONNER

This donning and doffing aid is a soft device that allows one to quickly and painlessly slide medical compression stockings over an arm or a leg. It is ideal for independently donning and doffing as well as for support through nursing staff. The cone simplifies the rolling of a compression stocking onto the DOFF N' DONNER.



Cone



DOFF N' DONNER (DnD)

ACCESSORIES

Care and additional products

The right care is essential to ensure that the stockings retain their pressure and shape for as long as possible and fit perfectly. Further additional products help to successfully master everyday life with compression stockings.

FOR BETTER HOLD AND DURABILITY OF GRIP TOPS: CLEAN PAD

Disposable cleaning wipe with conditioning ingredients for improved adhesion of the silicone coating and increased durability of the grip tops.

Does not contain alcohol.

1 box with 10 pads



FOR EASY WASHING OF THE STOCKINGS: WASHING SOLUTION SWS

Special hand wash for medical compression stockings and tights. The WASHING SOLUTION meets all the requirements for optimum stocking care and can also be used for all delicate laundry.

Daily washing allows the compression stockings to return to their original shape and compression state after wearing.


1 single bottle (250 ml) / 6 bottles à 250 ml



FOR SECURING STOCKINGS WITHOUT A GRIP TOP: FIX

A special roll-on adhesive for securing compression stockings to arms and legs without a grip top. It is gentle to the skin and follows all body movements. It can also be used to secure orthopedic and surgical aids etc. to the body. FIX can be washed off with water and leaves no stains.

1 piece (60 ml) / 12 pieces à 60 ml





"I want a straightforward compression solution that works. A large choice of colors is not important to me. That is why I choose Dynaven."



"My segment is Dynaven. A compression stocking must be practical and sturdy. No more and no less."

DYNAVEN

Dynaven products are functional and effective – no more and no less. They are suitable for pragmatic users who want a compression stocking that offers simple, solid and practical efficacy. Dynaven is our second product brand after Sigvaris.

Dynaven CLASSIC p. 47






Dynaven CLASSIC

For her and him

Dynaven CLASSIC is the functional, effective and opaque stocking for perfectly fitting basic care. It is ideal for both women and men. The Dynaven CLASSIC range is now also available as made-to-measure.

Series and made-to-measure range:

Dynaven CLASSIC		A-D Calf		A-G Thigh	A-T Pantyhose	A-T Maternity
Matchcode DCL		Soft top	Grip top	Grip top	Soft top	Soft top
CCL	1	●	●	●	●	●
	2	●	●	●	●	●
Toe		●	●	●	●	●
		●	●	●	●	●
Grip top	Quadrinova knobbed grip top		●	●		
Colors	 black	●	●	●	●	●
	 caramel	●	●	●	●	●

Complete overview by compression levels and indications / needs at the front inside the cover; information on special models / needs on **p.50**



TAKING MEASUREMENTS CORRECTLY

Among other things, medical compression stockings are characterized by having to measure the legs to determine the size. This is to ensure that the stocking fits, is comfortable to wear and that compression on the leg is correct.


When

Take measurements first in the morning whenever possible. The leg should be free of edema/swelling. Take measurements on the skin (not over trousers).

How

Circumference: Place the measuring tape firmly around the leg / hip / waist without constricting.

Length: Measure length at the inside of the leg, with leg straight (with the exception of lT).



Note:

The length measurement is taken from the floor to the corresponding measuring point (e.g. lD , lG , lT).

lA open toe
 lZ closed toe



A video on „Taking measurements correctly“ can be found here:

TIPS AND TRICKS

The following points are important when advising your customers:

Tip 1: Explaining the benefits of compression stockings

Show the customer exactly how important, effective and above all comfortable compression therapy can be. This is the only way to overcome existing prejudices and achieve the best possible therapy success.

Tip 2: Training donning and doffing techniques

Show the customer how to put on and take off the compression stocking correctly. Then allow the customer to put on and take off the stocking again on his/her own so that you can make any corrections if necessary.

Tip 3: Give instructions on handling and care

For your customers to feel completely comfortable with their compression stockings, it is important to give them some basic tips on handling and care.



From early till late:

The compression stockings should be worn from morning to evening or as prescribed by the doctor. The stockings must be put on after getting up before the legs start to swell.



Gloves help:

Our special gloves should always be worn when putting on and taking off compression stockings. This makes donning and doffing much easier and protects the stocking.



Donning and doffing aids as support:

Further support can be provided by alternative donning and doffing aids such as the sim-slide, the magnide on/off, the DOFF N' DONNER or MELANY.



Daily washing:

The stockings should be washed daily by hand or in the washing machine using a gentle cycle with a maximum temperature of 40 °C and a laundry net. The detergent should not contain a fabric softener.



Drying the stockings:

The stocking must not be dried in a tumble drier or on a heater. The clothesline is best suited for drying.

SPECIAL MODELS / SPECIAL REQUIREMENTS

As every customer has different requirements, we also offer optimal products for special needs in addition to our standard models, borders and sizes. It is important to respond individually to customers and to find the best possible product. Only a product that fits and in which one feels comfortable will be worn.

SPECIAL BORDERS

Our calf socks close with a soft top and our thigh stockings with a grip top. However, there are situations which require a different top.

Calf with grip top

The alternative calf model offers an even better grip on the leg should the calf slip with the soft top. The calf with grip top is mainly used when the calf section is only slightly developed (i.e. when the circumferential mass in the cC and cD range is similar).

Available as Style **SEMITRANSSPARENT**, Essential **COMFORTABLE, THERMOREGULATING & CLASSICAL**, Specialities **TRADITIONAL** and Dynaven **CLASSIC**.

Thigh with soft top

This thigh stocking variant with soft top is used for customers with sensitive or allergy-prone skin. The special adhesive FIX can be used to provide the stockings with optimum hold.

Available as Essential **COMFORTABLE & THERMOREGULATING** and Specialities **TRADITIONAL**.

Thigh with waist attachment

This model variant is preferred by customers who do not want to wear a grip top and enables improved quality of life in old age (no tedious donning and doffing in case of age-related complaints such as incontinence etc.). Can also be worn as a left/right combination instead of a pantyhose. In the waist area with a strap with adhesive fastener.

Available as Essential **COMFORTABLE, THERMOREGULATING & CLASSICAL** and Specialities **TRADITIONAL**.

SPECIAL SIZES

Ensuring that the stocking fits perfectly is extremely important for compliance and the effectiveness of the stockings. In addition to our standard sizes, which already cover a wide range of sizes, we can offer individual solutions with short lengths or special stockings for customers with large feet.

Short lengths

This model variant is suitable for small persons or persons with a shorter lower leg or leg length.

Short length ℓD Calf: 31–34 cm

Short length ℓG Pantyhose: 56–64 cm

Available as Essential **THERMOREGULATING**.

Large feet

Here we distinguish between voluminous feet and feet with excessive length.

- Voluminous foot:
The Active **MASCULINE** stocking is perfect for large feet in all sizes. It is specifically designed for the male foot and can also relieve symptoms associated with hallux valgus (deformity of the big toe). Suitable for large ball circumferences (cA), extensive heel instep.
- Foot with overlength:
For the Essential **THERMOREGULATING** maxifoot product, we offer a medical compression stocking with an enlarged foot section.

SPECIAL MODELS

In addition to our standard models, we also offer various other models to meet as many needs as possible.

Half-thigh stockings

This thigh variant is ideal for customers for whom varicose veins in the back of the knee make a calf sock unsuitable and who do not wish to wear a pantyhose or conventional thigh stocking.

Available as Specialities **TRADITIONAL** and Essential **CLASSICAL**.

Pantyhose Bodyform

With the pantyhose Bodyform compression extends into the body region. In normal compression pantyhoses, the compression pressure ends at the level of the cG measurement. The Bodyform pantyhose is therefore suitable for persons with soft connective tissue and supports the silhouette of the hip. A further medical indication for compression in the body region is lipedema in the hip area.

Available as Essential **COMFORTABLE**.



Examples of special tops

TOE OPTIONS

Most of our compression stockings come with a choice of two toe options. In most cases the choice can be left to the customer.

Open foot

- For customers who do not require a closed toe.
- More pleasant to wear at warm temperatures.
- Easier to put on, therefore ideal for elderly people.
- Allows checking for possible discoloration and pressure points of the toes.



Toe options

Closed toe

- Required for medical reasons if compression pressure is required in the toe area, e.g. for edemas in the foot area.
- For orthopedic reasons with a strongly pronounced hallux.

COMPRESSION VALUES, INDICATIONS AND CONTRAINDICATIONS

Wearing compression stockings is an effective therapy against chronic and acute venous diseases, for preventing travel thrombosis or spider veins or varicose veins during pregnancy. We recommend using medical compression stockings only after consultation with a (specialist) doctor.

COMPRESSION VALUES IN THE ANKLE AREA

There are four compression classes. Experts also speak of CCL 1 to 4: the higher the class, the stronger the compression: this is the pressure exerted by the medical compression stocking. Pressure is expressed in either mmHg (millimeters of mercury) or kPa (kilopascals).

- **Compression level 1 / CCL 1:**
18–21 mmHg; kPa value: 2.4–2.8
- **Compression level 2 / CCL 2:**
23–32 mmHg; kPa value: 3.1–4.3
- **Compression level 3 / CCL 3:**
34–46 mmHg; kPa value: 4.5–6.1
- **Compression level 4 / CCL 4:**
min. 49 mmHg; kPa value: ≥6.5

INDICATIONS

Compression therapy is generally prescribed by a doctor for the following diseases¹:

Chronic venous diseases

- Improving venous symptoms
- Improving quality of life for chronic venous diseases
- Preventing and treating venous edema
- Preventing and treating venous skin changes
- Eczema and pigmentation
- Lipodermatosclerosis and atrophie blanche
- Treating venous leg ulcers
- Treating mixed (arterial and venous) leg ulcers (taking the contraindications into account)
- Preventing recurrent venous leg ulcers
- Pain reduction with venous leg ulcers
- Varicosis
- Initial phase after varicose therapy
- Functional venous insufficiency (with obesity, and seated and standing occupations)
- Venous malformations

Venous thrombotic diseases

- Superficial venous thrombosis
- Deep vein thrombosis
- Arm vein thrombosis
- Post-thrombotic conditions
- Post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients

Edema

- Lymphedema
- Edema during pregnancy
- Post-traumatic edema
- Post-operative edema
- Post-operative reperfusion edema
- Idiopathic cyclic edema
- Lipedema from stage II
- Congestion as a result of immobility (arthrogenic congestion syndrome, paresis and partial paresis of the extremities)
- Occupational edema (standing and seated occupations)
- Drug-induced edema, if switching drugs is not possible

Other indications

- Obesity with functional venous insufficiency
- Inflammatory dermatosis of the legs
- Nausea and dizziness during pregnancy
- Congestive conditions during pregnancy
- Post-burn conditions
- Scar treatment

CONTRAINDICATIONS

Medical compression stockings must not be worn in the following cases:

- Advanced peripheral arterial occlusive disease (if one of the following parameters applies: ABI <0.5, ankle blood pressure <60 mmHg, toe blood pressure <30 mmHg or TcPO₂ <20 mmHg instep).
- Decompensated cardiac insufficiency (NYHA III + IV)
- Septic thrombophlebitis
- Phlegmasia cerulea dolens

Relative contraindications, whereby benefits and risks must be weighed up, in particular:

- Severe exuding dermatosis
- Intolerance to compression material
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g., in diabetes mellitus)
- Primary chronic polyarthritis

RISKS AND SIDE EFFECTS

MCS can cause skin necrosis and pressure damage to the peripheral nerves especially if handled improperly. In cases of sensitive skin, itching, scaling, and signs of inflammation may occur under the compression aids. This is why it is advisable to use adequate skin care under the compression aids.




¹ Guideline of the German Society of Phlebology on Medical Compression Stockings Phlebology 2019; AWMF (Association of the Scientific Medical Societies) registration: 037/005


SIZING SYSTEMS

Circumference and length size system for Essential **COMFORTABLE**,
Essential **THERMOREGULATING**, Style **SEMITRANSSPARENT**, Active **MASCULINE**

		Circumference cm											
		XXS*	XXS Plus*	XS	XS Plus	S	S Plus	M	M Plus	L	L Plus	XL	XL Plus
cT		60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130
cH		75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
cG	without grip top	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	73-85
	with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
cF		38-49	41-52	41-52	44-55	44-55	49-60	49-60	53-64	53-64	59-70	59-70	63-78
cE		29-36	33-40	33-40	37-44	37-44	39-46	39-46	41-48	41-48	45-52	45-52	46-53
cD		23-30	27-34	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48	41-48	42-50
cC		25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
cB₁		22-26	22-26	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37	33-39	33-39
cB		18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29	29-32	29-32
cY		27-31	27-31	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38	33-39	33-39
cA		18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29	29-32	29-32



Length cm					
short*	31-34	short	56-64	short**	56-64
normal	35-40	normal	65-74	normal	65-74
long	41-45	long	75-84	long	75-84




* short length only available in Essential THERMOREGULATING A-D with soft top in CCL 2, colors black and nature; not available as maxifoot
** short length only available in Essential THERMOREGULATING A-T with soft top in CCL2 und 3, color nature


A-D = Calf; **A-G** = Thigh; **A-T** = Pantyhose; **A-T** Maternity = for pregnant and oversized women; circumferential measurements cH and cT are replaced by a special adapting trouser component

Circumference and length size system for Essential **CLASSICAL**

		Circumference cm					
		S	S Plus	M	M Plus	L	L Plus
cT		up to 80	up to 80	up to 90	up to 90	up to 100	up to 100
cH		up to 110	up to 110	up to 120	up to 120	up to 130	up to 130
cG		42-52	50-60	50-60	58-68	58-68	66-76
cF		41-50	45-54	48-57	53-62	55-64	61-70
cE		32-37	36-41	35-40	40-45	38-43	44-49
cD		28-34	32-38	31-37	36-42	34-40	40-46
cC		31-36	35-40	34-39	39-44	37-42	43-48
cB₁		24-29	24-29	29-32	29-32	32-36	32-36
cB		19-23	19-23	23-26	23-26	26-29	26-29
cY		28-32	28-32	32-35	32-35	35-38	35-38
cA		19-23	19-23	23-26	23-26	26-29	26-29



Length cm							
short	33-37	short	56-60	short	61-70	short	61-70
long	38-44	long	61-70	long	71-80	long	71-80




A-D = Calf; **A-F** = Half-thigh; **A-G** = Thigh; **A-T** = Pantyhose


SIZING SYSTEMS

Circumference and length size system for Specialities **TRADITIONAL**, Specialities **ULCER X**

	Circumference cm								
	XS	S	S Plus	M	M Plus	L	L Plus	XL	XL Plus
cT		up to 80	up to 80	up to 90	up to 90	up to 100	up to 100		
cH		up to 110	up to 110	up to 120	up to 120	up to 130	up to 130		
cG	40-50	42-52	50-60	50-60	58-68	58-68	66-76	61-71	69-79
cF	39-48	41-50	45-54	48-57	53-62	55-64	61-70	58-67	64-73
cE	29-34	32-37	36-41	35-40	40-45	38-43	44-49	41-46	48-53
cD	26-32	28-34	32-38	31-37	36-42	34-40	40-46	38-43	44-49
cC	28-31	31-36	35-40	34-39	39-44	37-42	43-48	40-45	47-52
cB₁	21-24	24-29	24-29	29-32	29-32	32-36	32-36	35-39	35-39
cB	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32
cY	25-28	28-32	28-32	32-35	32-35	35-38	35-38	38-41	38-41
cA	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32



Length cm							
short	33-37	short	56-60	short	61-70	short	61-70
long	38-44	long	61-70	long	71-80	long	71-80
				x-long*	81-90		



* x-long length only available for TRADITIONAL A-G Post-OP

A-D = Calf; **A-F** = Half-thigh; **A-G** = Thigh; **A-T** = Pantyhose; **A-T Maternity** = for pregnant and oversized women: circumferential measurements cH and cT are replaced by a special adapting trouser component

Available sizes for Specialities **TRADITIONAL**, Specialities **ULCER X**

	Compression class	Sizes									Length			
		XS	S	S PLUS	M	M PLUS	L	L PLUS	XL	XL PLUS	short	long	x-long	
A-D Calf	CCL 2	•	•	•	•	•	•	•	•	•	•	•	•	
	CCL 3		•	•	•	•	•	•	•	•	•	•	•	
	CCL 4		• ¹	• ¹	• ¹	• ¹	• ¹	• ¹	• ¹	• ¹	•	•	•	
A-F Half-thigh	CCL 2	• ^{1,2}	•	•	•	•	•	•	•	• ^{1,2}	• ^{1,2}	•	•	
	CCL 3		•	•	•	•	•	•	•	•	•	•	•	
	CCL 4			• ¹		• ¹		• ¹		• ¹		•	•	
A-G Thigh	CCL 2	• ^{1,2}	•	•	•	•	•	•	•	• ^{1,2}	• ^{1,2}	•	•	
	CCL 3		•	•	•	•	•	•	•			•	•	
	CCL 4			• ¹		• ¹		• ¹		• ¹		•	•	
A-G Thigh with waist attachment	CCL 2		•	•	•	•	•	•	•			•	•	
	CCL 3			•		•		•		•		•	•	
A-G Post-OP	CCL 2			•		•		•				•	•	•
	CCL 3			•		•		•				•	•	•
A-T Pantyhose A-T Maternity	CCL 2			•		•		•				•	•	
A-T for men	CCL 2			•		•		•				•	•	
ULCER X Kit*	~45mmHg (~6.0 kPa)	•	•	•	•	•	•	•	•	•	•	•	•	

¹ model only available without grip top
² model only available in long

* Notes on sizes of the separately available understockings for **ULCER X** specialities:


Size ULCER X Kit	Size understocking (set of 4)
XS short	S short
XS long	S long
S short, S Plus short	S short
S long, S Plus long	S long
M short, M Plus short	M short
M long, M Plus long	M long
L short, L Plus short	L short
L long, L Plus long	L long
XL short, XL Plus short	XL short
XL long, XL Plus long	XL long

SIZING SYSTEMS

Circumference and length size system for Specialities **DIABETIC COMPRESSION SOCKS**

	Circumference cm			
	S	M	L	XL*
cD	27-33	31-39	35-45	41-51
cB	18-22	22-26	26-29	29-36

* Size XL only available for men




Length cm	
short	≤ 41
long	> 41



Shoe sizes (guide values)			
		short	long
Euro	women	up to 39	up to 42
	men	up to 43	up to 46
UK	women	up to 6	up to 8
	men	up to 9	up to 11
US	women	up to 8.5	up to 10.5
	men	up to 10	up to 12
AUS/NZL	women	up to 7	up to 8.5
	men	up to 8.5	up to 11.5

The two variants DIABETIC COMPRESSION SOCKS women und DIABETIC COMPRESSION SOCKS men differ in the size of the foot area. In addition to the standard sizing, the guide values can be consulted for exceptionally large or small shoe sizes. In case of doubt, the patient can change alternatively to men's or to women's socks

Circumference and length size systems for Specialities **THROMBO & THROMBO X**



THROMBO X Circumference (cm)	cB	cD	cG
small	18-22 cm	29-35 cm	43-54 cm
medium	22-26 cm	34-40 cm	48-60 cm
large	26-29 cm	39-45 cm	54-66 cm
x-large	29-32 cm	44-51 cm	60-72 cm
xx-large ¹	32-35 cm	50-56 cm	67-80 cm

¹ only for length normal


THROMBO X Length (cm)	lD	lG
short	33-38 cm	60-70 cm
normal	38-43 cm	70-80 cm
long	43-48 cm	80-90 cm

THROMBO Circumference (cm)	cB	cD	cG
x-small ²	17-19 cm	≤ 34 cm	≤ 52 cm
small	19-23 cm	≤ 38 cm	≤ 60 cm
medium	23-26 cm	≤ 42 cm	≤ 68 cm
large	26-29 cm	≤ 46 cm	≤ 76 cm

² not available for A-G Thigh with waist attachment

THROMBO Length (cm)	lD	lG
short	≤ 37 cm	≤ 70 cm
long	≥ 38 cm	≥ 71 cm

Circumference and length size system for Specialities **ARM SLEEVES**



¹ only for mitten (A-C)
² no Plus sizes available for ADVANCE A-G with soft top and C-G with soft top

	ADVANCE Circumference cm										TRADITIONAL Circumference cm		
	XS	XS Plus ²	S	S Plus ²	M	M Plus ²	L	L Plus ²	XL	XL Plus ²	S	M	L
cF	24-32	32-38	26-34	34-40	29-37	37-43	32-40	40-46	34-42	42-48	24-30	30-36	36-42
cD											20-24	24-28	28-32
cC	13-14	13-14	14-16	14-16	16-18	16-18	18-20	18-20	20-23	20-23	13-16	16-19	19-22
cA¹	16-20	16-20	17-21	17-21	18-22	18-22	19-23	19-23	20-24	20-24	16-19	19-22	22-25

ADVANCE Length cm	
normal	< 43
long	> 43

A-G	A-G	C-G	C-G
with mitten with soft top	with mitten with Sensinova grip top	without mitten with soft top	without mitten with Sensinova grip top


TRADITIONAL Length cm	
normal	40 - 46
long	46 - 52

A-C	C-H	C-G	A-H	A-G
separate mitten	without mitten with support straps (arm)	without mitten with grip top	with mitten with support straps (arm)	with mitten with grip top

SIZING SYSTEMS

Circumference and length size system for Dynaven **CLASSIC**

	Circumference cm					
	S	S Plus	M	M Plus	L	L Plus
cT	up to 80	up to 80	up to 90	up to 90	up to 100	up to 100
cH	up to 110	up to 110	up to 120	up to 120	up to 130	up to 130
cG	42-52	50-60	50-60	58-68	58-68	66-76
cF	41-50	45-54	48-57	53-62	55-64	61-70
cE	32-37	36-41	35-40	40-45	38-43	44-49
cD	28-34	32-38	31-37	36-42	34-40	40-46
cC	31-36	35-40	34-39	39-44	37-42	43-48
cB₁	24-29	24-29	29-32	29-32	32-36	32-36
cB	19-23	19-23	23-26	23-26	26-29	26-29
cY	28-32	28-32	32-35	32-35	35-38	35-38
cA	19-23	19-23	23-26	23-26	26-29	26-29



Lenght cm					
short	33 – 37	short	61 – 70	short	61 – 70
long	38 – 44	long	71 – 80	long	71 – 80

A-D = Calf; **A-G** = Thigh; **A-T** = Pantyhose; **A-T Maternity** = for pregnant and oversized women: circumferential measurements **cH** and **cT** are replaced by a special adapting trouser component

WE ARE DEDICATED TO HELPING PEOPLE
FEEL THEIR BEST. EVERY DAY.

Export

SIGVARIS AG
Gröblistrasse 8
9014 St. Gallen
Switzerland
Phone +41 71 272 40 00
Fax +41 71 272 40 01

Great Britain

SIGVARIS Britain LTD.
1 Imperial Court
Magellan Close
Andover
SP10 5NT
Phone +44 1264 326 666
Fax +44 1264 369 950

Australia / New Zealand

SIGVARIS (Australia) Pty. Ltd.
8/53 Stanley Street
West Melbourne VIC 3003
Phone 1300 SIGVARIS (1300 744 827)
Fax +61 39329 3773
info.au@sigvaris.com

www.sigvaris.com

printed in
switzerland

