

sigvaris

Essential THERMOREGULATING



sensitive skin

Soft: is also suitable for sensitive legs due to 14 % cotton content and soft top



thermoregulating

Cooling: evaporation of moisture stored in the cotton, resulting in a cooling effect










Warming: the air in the spaces between the knitted fabric and the fibers forms an insulating layer at cold temperatures



opaque

Attractive: discreet, non-transparent ribbed structure, conceals venous problem zones perfectly



Matchcode ETH		 A-D Calf		 A-G Thigh			 A-T Pantyhose	 A-T Maternity
		Soft top ^{1,2}	Grip top	Soft top	Grip top	Waist attachment	Soft top ³	Soft top
CCL	1	●		●	●	●	●	●
	2	●	●	●	●	●	●	●
	3	●			●	●	●	
Toe		●	●	●	●	●	●	●
		●	●	●	●	●	●	●
Grip top	Sensinova grip top				● CCL 1/2			
	Quadrinova knobbed grip top		●		●			
	Micro knobbed grip top				● CCL 1/2			
Colors	 black	●	●	●	●		●	●
	 nature	●	●	●	●	●	●	●
	 marine	available as A-D with soft top, CCL 2, closed toe, normal and long						


¹ CCL 2 with closed toe is also available as maxifoot in black or nature, without XXS and XXS Plus

² CCL 2 with open and closed toe, is also available as short length in black or nature


³ CCL 2 and 3 with open and closed toe, is also available as short length in nature



		Circumference cm											
		XXS	XXS Plus	XS	XS Plus	S	S Plus	M	M Plus	L	L Plus	XL	XL Plus
cT		60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130
cH		75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
cG	without grip top	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	73-85
	with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
cF		38-49	41-52	41-52	44-55	44-55	49-60	49-60	53-64	53-64	59-70	59-70	63-78
cE		29-36	33-40	33-40	37-44	37-44	39-46	39-46	41-48	41-48	45-52	45-52	46-53
cD		23-30	27-34	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48	41-48	42-50
cC		25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
cB₁		22-26	22-26	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37	33-39	33-39
cB		18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29	29-32	29-32
cY		27-31	27-31	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38	33-39	33-39
cA		18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29	29-32	29-32



Length cm					
short*	31 - 34	short	56 - 64	short**	56 - 64
normal	35 - 40	normal	65 - 74	normal	65 - 74
long	41 - 45	long	75 - 84	long	75 - 84



* short length only available in Essential THERMOREGULATING A-D with soft top in CCL 2, colors black and nature; not available as maxifoot

** short length only available in Essential THERMOREGULATING A-T with soft top in CCL 2 and 3, color nature

A-D = Calf; **A-G** = Thigh; **A-T** = Pantyhose; **A-T** Maternity = for pregnant and oversized women: circumferential measurements cH and cT are replaced by a special adapting trouser component

Export • SIGVARIS AG • Gröblistrasse 8 • 9014 St. Gallen, Switzerland • Phone +41 71 272 40 00 • Fax +41 71 272 40 01

Great Britain • SIGVARIS Britain LTD. • 1 Imperial Court • Magellan Close • Andover, SP10 5NT • Phone +44 1264 326 666 • Fax +44 1264 369 950

Australia / New Zealand • SIGVARIS (Australia) Pty. Ltd. • 8/53 Stanley Street • West Melbourne VIC 3003 • Phone 1300 SIGVARIS (1300 744 827) • Fax +61 39329 3773 • info.au@sigvaris.com

www.sigvaris.com