

sigvaris



TRADITIONAL

Maximum patient acceptance
and optimal therapy results

Live your life. Be yourself.

This motto guides us.

It is aimed at the wearers of our compression solutions. They are our focus. Individuals with diverse medical and personal needs. People who only want one thing - to lead an active life while remaining true to their style. Without compromises.

We want to develop the best compression solutions for people. Each and every individual should feel good. Every day. We have been focusing on this **for over 150 years.**

Our Premium portfolio of medical compression stockings is made up of the following **segments** that address various needs and indications:

- **Essential**
- **Style**
- **Active**
- **Specialities**

Each segment is designed to not only let you find a suitable, but also the best compression solution. On the one hand, the focus is on the effect and medical aspects, on the other, the products should make the customer feel good and help to lead a life without restrictions.

Specialities

Every person is unique and their requirements for compression solutions can change depending on their life situation.

However, sometimes the medical aspect is of primary concern and the needs of the stocking wearer cannot be met with standard compression solutions. These are situations where products from the Specialities segment come into their own. They are characterized in that **each stocking has been developed for a very specific medical indication or special need.**

In this folder, we take an in-depth look at **TRADITIONAL**, which can be found in the Specialities segment.

TRADITIONAL is manufactured in St. Gallen and therefore stands for **Quality made in Switzerland.** This is something we are proud of.

Specialities

TRADITIONAL



For her and him

TRADITIONAL is the number one choice for when medical efficacy becomes a top priority.



The yarn's unique stress/strain behavior guarantees excellent pressure stability throughout the day and the effective treatment of severe venous diseases. This medical compression stocking is also easy to slip on and off in high pressure classes thanks to its quite stretchable materials. It is robust, durable and breathable.

Recommended by physicians worldwide

Medical benefits confirmed by numerous scientific studies (www.stemmerlibrary.com).




Product range overview TRADITIONAL

TRADITIONAL Matchcode TRA		A-D Calf		A-F Half-thigh		A-G Thigh		
		Soft top	Grip top	Soft top	Grip top	Soft top	Grip top	Waist attachment
CCL	2	•	•	•	•	•	•	•
	3	•	•	•	•	•	•	•
	4	•		•		•		
Toe		•	•	•	•	•	•	•
Grip top	Knobbed grip top		•		•		•	
Colors	 beige	•	•	•	•	•	•	•

* Kit compromises: 2x understocking (without compression), 1x TRADITIONAL

Product features




Durable material

High-quality materials result in less wear and tear. TRADITIONAL is therefore ideal for demanding everyday routines. Wear tests have proven: TRADITIONAL has up to 6 times better* abrasion results in direct comparison with jeans (* = less wear).



Knobbed grip top

The round smooth knobbes are particularly pleasant on the skin and do not leave any marks. With soft protective seam edge.




High breathability

The wide-meshed structure of the knitted fabric ensures high air permeability.



Soft top

A particularly soft top for customers with sensitive or allergy-prone skin.







Natural rubber

Best possible elastic recovery plus very high pressure stability.

Choice of colors



 A-G Post-OP	 A-T Strumpfhose	 A-T Maternity	 A-T For men
Waist attachment*	Soft top	Soft top	Soft top
•	•	•	•
•			
•	•	•	•
•	•	•	•

Special models

Half-high stockings

This thigh variant is ideal for customers for whom varicose veins in the back of the knee make a calf sock unsuitable and who do not wish to wear a pantyhose or conventional thigh stocking.

Post surgery

This model is suitable for patients who have undergone leg surgery. Kit comprises:
 2x understocking (without compression)
 1x TRADITIONAL A-G with waist attachment


Additional understockings are available separately.

More information about the other special models can be found in the assortment brochure.

Sizing systems

Circumference and length for **Specialities TRADITIONAL**

	Circumference cm								
	XS	S	S Plus	M	M Plus	L	L Plus	XL	XL Plus
cT		up to 80	up to 80	up to 90	up to 90	up to 100	up to 100		
cH		up to 110	up to 110	up to 120	up to 120	up to 130	up to 130		
cG	40-50	42-52	50-60	50-60	58-68	58-68	66-76	61-71	69-79
cF	39-48	41-50	45-54	48-57	53-62	55-64	61-70	58-67	64-73
cE	29-34	32-37	36-41	35-40	40-45	38-43	44-49	41-46	48-53
cD	26-32	28-34	32-38	31-37	36-42	34-40	40-46	38-43	44-49
cC	28-31	31-36	35-40	34-39	39-44	37-42	43-48	40-45	47-52
cB ₁	21-24	24-29	24-29	29-32	29-32	32-36	32-36	35-39	35-39
cB	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32
cY	25-28	28-32	28-32	32-35	32-35	35-38	35-38	38-41	38-41
cA	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32










Length cm							
short	33 - 37	short	56 - 60	short	61 - 70	short	61 - 70
long	38 - 44	long	61 - 70	long	71 - 80	long	71 - 80
				x-long*	81 - 90		

* x-long length only available for TRADITIONAL A-G Post-OP

A-D = Calf; **A-F** = Half-thigh; **A-G** = Thigh; **A-T** = Pantyhose; **A-T** Maternity = for pregnant and oversized women: circumferential measurements cH and cT are replaced by a special adapting trouser component

Should you not find what you are looking for in our extensive TRADITIONAL standard range, we would be pleased to supply you with our custom-made TRADITIONAL products.

Available sizes for Specialties TRADITIONAL

	Compression class	Sizes									Length			
		XS	S	S PLUS	M	M PLUS	L	L PLUS	XL	XL PLUS	short	long	x-long	
 A-D Calf	CCL 2	●	●	●	●	●	●	●	●	●	●	●	●	
	CCL 3		●	●	●	●	●	●	●	●	●	●	●	
	CCL 4		● ¹	● ¹	● ¹	● ¹	● ¹	● ¹				●	●	
 A-F Half-thigh	CCL 2	● ¹	●	●	●	●	●	●	●	● ¹	● ¹	●	●	
	CCL 3		●	●	●	●	●	●				●	●	
	CCL 4			● ¹		● ¹		● ¹				●	●	
 A-G Thigh	CCL 2	● ^{1,2}	●	●	●	●	●	●	●	● ^{1,2}	● ^{1,2}	●	●	
	CCL 3		●	●	●	●	●	●				●	●	
	CCL 4			● ¹		● ¹		● ¹					●	
 A-G Thigh with waist attachment	CCL 2		●	●	●	●	●	●				●	●	
	CCL 3			●		●		●				●	●	
 A-G Post-OP	CCL 2			●		●		●				●	●	●
	CCL 3			●		●		●				●	●	●
 A-T Pantyhose A-T Maternity	CCL 2			●		●		●				●	●	
 A-T for men	CCL 2			●		●		●				●	●	

¹ model only available without grip top, ² model only available in long

Care

The stockings should be washed by hand after wearing or at max. 40 °C with the gentle wash cycle in the washing machine in a laundry net. It is important to use a mild detergent without fabric softener. The Sigvaris Washing Solution is suitable for washing by hand.

To preserve maximum durability of the product, please inform your customers that we do not recommend using a dryer/heater. The stockings must not be ironed or dry cleaned. The stockings should be stored dry at room temperature away from direct sunlight.



TRADITIONAL contains latex.

Donning and doffing

Donning and doffing of compression stockings presents a challenge for some wearers.



Recommend the **practical donning and doffing videos** on our website to your customers.

These labels commit us to the highest quality

We fulfil all quality and test regulations of the RAL Deutsches Institut für Gütesicherung und Kennzeichnung e.V. [German Institute for Quality Assurance and Labeling] for medical compression stockings. Our medical compression stockings have been tested for human ecology and comply with the internationally binding Eco-Tex standard 100. Our medical compression stockings meet the basic safety and performance requirements of the relevant EU regulations.

sigvaris

Export

SIGVARIS AG
Gröblistrasse 8
9014 St. Gallen
Switzerland
Phone +41 71 272 40 00
Fax +41 71 272 40 01

Great Britain

SIGVARIS Britain LTD.
1 Imperial Court
Magellan Close
Andover
SP10 5NT
Phone +44 1264 326 666
Fax +44 1264 369 950

Australia / New Zealand

SIGVARIS (Australia) Pty. Ltd.
8/53 Stanley Street
West Melbourne VIC 3003
Phone 1300 SIGVARIS (1300 744 827)
Fax +61 39329 3773
info.au@sigvaris.com

www.sigvaris.com

**We are dedicated to helping people
feel their best. Every day.**